

Disinfectant Spray Q (Lemon Scent) - 17oz. Bottle

Product Spec Sheet



Provided by

MyBinding.com
When Image Matters.

Call Us at 1-800-944-4573

CL1002 TECHNICAL DATA SHEET

CLAIRE DISINFECTANT SPRAY Q LEMON SCENT



Claire Disinfectant Spray Q EPA Reg. No. 706-111 has demonstrated effectiveness against viruses similar to 2019 Novel Coronavirus (2019-nCoV) on hard non-porous surfaces. Therefore, Claire Disinfectant Spray Q can be used against 2019 Novel Coronavirus (2019-nCoV) when used in accordance with the directions for use against Poliovirus type 1 on hard, non-porous surfaces.

Refer to the CDC website <https://www.cdc.gov/coronavirus/2019-ncov/index.html> for additional information.

PRODUCT DESCRIPTION:

- Hospital grade multipurpose disinfectant for hard, non-porous surfaces
- Soft spot surface sanitizer
- Controls mold and mildew
- Eliminates odors
- Kills 99.9% germs
- Bactericidal, Antibacterial, Fungicidal, *Virucidal, Tuberculocidal

Disinfects and kills the following bacteria and *Viruses in 3 minutes on pre-cleaned hard, non-porous surfaces:

Kills Bacteria: *Pseudomonas aeruginosa* (ATCC 15442), *Salmonella enterica* (ATCC 10708), *Staphylococcus aureus* (ATCC 6538), *Acinetobacter baumannii* Multi Drug Resistant (ATCC 19606), Carbapenem Resistant *Escherichia coli* (CDC 81371), Carbapenem Resistant *Klebsiella pneumoniae* (ATCC BAA-1705), *Enterobacter aerogenes* (ATCC 13048), *Enterococcus faecium* (ATCC 51559), *Escherichia coli* (ATCC 11229), *Klebsiella pneumoniae* (ATCC 4352), *Listeria monocytogenes* (ATCC 19117), *Staphylococcus aureus*- MRSA (ATCC 33592), *Staphylococcus aureus*- CA-MRSA Genotype USA 300 (NARSA NRS 384), *Staphylococcus aureus*- CA-MRSA Genotype USA 400 (NARSA NRS 123), *Shigella dysenteriae* (ATCC 11835), *Streptococcus pyogenes* (ATCC 19615), *Staphylococcus aureus* -VISA (CDC HIP 5836), *Enterococcus faecalis* -VRE (ATCC 51575), *Campylobacter jejuni* (ATCC 29428), *Candida albicans* (ATCC 10231), *Mycobacterium bovis* (TB) (10 minute contact time).

***Viruses:** Adenovirus Type 2*, 2009-H1N1 Influenza A virus*, Influenza B virus*, Human Immunodeficiency virus type 1 (HIV-1)*, Herpes simplex virus type 1*, Herpes simplex virus type 2*, Rotavirus*, Vaccinia virus*, Rhinovirus type 37*.

Poliovirus type 1 in 5 minutes on pre-cleaned hard, non-porous surfaces.

Fungi: *Trichophyton interdigitale* (Athlete's Foot Fungus) (ATCC 9533) (10 minute contact time).

USE IN:

Healthcare, schools, hotels, restaurants, commercial establishments. Use in homes on toilets, sinks, showers, bathtubs, doorknobs, light switches, garbage cans, and high touch areas.

Sanitizes soft non-food contact surfaces such as upholstered furniture, stuffed pillows, draperies, etc.

SPECS:

Can Size: 20 oz.
Net Weight: 17 oz.
Shipping Weight: 17 Lbs.
Packaged: 12 cans per case
UPC Number: 7 13014 11002 2
GTIN: 10713014610024
Product Number: CL1002
Extender Tube: N

REGULATORY:

Ozone Depleting Compounds: None
Recyclable Package: Yes
VOC Compliant CARB: Yes
VOC Compliant OTC: Yes
NSF Certified: Yes
Flammability: Not required to be labeled as flammable

PHYSICAL CHARACTERISTICS:

Color: Colorless
Fragrance: Lemon
Shelf Life: 1 year +

See label side panel for First Aid and additional Precautionary Statements

24 HOUR MEDICAL EMERGENCY NUMBER:
1-866-836-8855

KEEP OUT OF REACH OF CHILDREN
CONSULT SDS BEFORE USING

